wellion[®]

MESH-NEBULIZER

User manual

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
MESH-INHALÁTOR Návod na použitie

Statements

- Thanks for purchasing the product.
- To ensure correct usage, please read the User Manual carefully before using this product.
- Please keep the User Manual handy for future reference.
- The company takes no responsibilities and does not provide free maintenance for any abnormal phenomena or damage due to incorrect usage, maintenance and / or storage.
- The company reserves final explanation rights to this manual.

1. Precautions

Please read the user manual carefully in order to ensure safe use.

Warning:

- Indicating a potentially hazardous situation that, if not avoided, could lead to serious consequences, including injury or death.

Attention:

- Emphasizing important information, instructions or explanations for better use.

Warnings:

- Pentamidine drugs, oily drugs and lipid-containing drugs are not applicable. In order to avoid damage to the nebulizing mesh, it is recommended to use standard nebulization drugs.

- Please follow doctors' advice regarding medication type, dosage and usage. Otherwise it may cause symptomatic deterioration.
- Do not use the device in the circuit of a respiratory anesthesia system and ventilator system (life-supporting system), otherwise it may cause incorrect operation.
- Please follow the specified operation instructions in this user manual, otherwise it may cause malfunctioning.
- Accessories of this device are not designed for reuse without prior cleaning, otherwise it may cause cross infection.
- When using this device for the first time or if the medication cup has not been used for a long time, both medication cup and mask must be cleaned and disinfected. Otherwise, it may cause bacterial infection.
- Do not use sodium hypochlorite disinfectant as this may lead to corrosion of the nebulizing mesh.
- To avoid damage do not use microwave ovens, tableware dryers or hair dryers to dry or disinfect the accessories.
- Please do not use water for nebulization, otherwise it may worsen your symptoms.
- For replacement please purchase only accessories with a registration certificate, otherwise it may cause allergic reactions and worsen symptoms.
- Please clean the accessories after disinfection, otherwise patient may inhale the residual disinfectant, which may cause symptomatic deterioration.
- Used medication must not be reused, please use fresh medication for every treatment. Otherwise patient may be infected by varieties of bacteria, causing symptomatic deterioration.
- Do not use the device to inhale water, otherwise it may cause symptomatic deterioration.

- Do not use the device at ambient temperature above 40 °C. Otherwise it may cause nasal mucosa injury or device failure.
- Do not rinse the main unit with water, do not immerse it in water and do not store the device in a humid environment. Otherwise it may cause device failure.
- Please clean the device each time after use, and dry it immediately after cleaning. Otherwise patient may be infected by varieties of bacteria.
- Please keep the device out of the reach of children and people with mental illness. Otherwise there is a risk of small parts being swallowed.
- Do not use the device near flammable or explosive gas or any anesthetic mixture. Otherwise it may cause personal injury.
- Be aware of possible suffocation danger if power cord wraps around children's neck.
- To ensure the device's safety, functionality and compliance to relevant regulations its unauthorized modification or use for any other purpose than the one defined is not allowed.

 Attention:

- If the device does not shutdown automatically when medication is used up, please immediately press the ON/OFF button to turn it off, in order to avoid damage to the nebulizing sheet. Refer to Chapter 6 Troubleshooting.
 - Clean the medication cup after each use. Otherwise, the device will not work properly.
 - Recommended maximum liquid capacity is 10ml, minimum capacity is 2ml.
 - During use temperature of the medication cup should not exceed 40 °C.
 - Do not rinse the medication cup under running water while it is still connected to the main body. This could result in water ingress.
 - Do not use this product near high-frequency electromagnetic transmitters and other high-frequency electronic products.
- 4

- Keep the device as upright as possible during use.
- Avoid dropping the device body or the medication cup. Do not subject them to severe shocks and blows.
- Do not touch the nebulizing mesh with a cotton swab or any sharp objects as this could lead to malfunctioning.
- This device should not be used by patients having facial/oral contusions, burns, inflammation or any kind of trauma. If any discomfort arises during use, please stop the treatment immediately and consult a doctor.
- Commonly used atomized drugs include moistening expectorants, bronchodilators and antibiotics, such as terbutaline sulphate solution for nebulization and ipratropium bromide solution for inhalation. The active ingredient must be in a liquid solution, either dissolved or suspended in a suitable solvent, allowing for effective nebulization. It must be compatible with the nebulizer device and the patient's respiratory system, avoiding irritation or adverse reactions. The solution's tonicity (salt concentration) and pH should be appropriate for safe and effective delivery to the respiratory tract, avoiding irritation or bronchoconstriction.
- Masks are made of plasticizer-free PVC material, medical professionals should refer to drug instruction leaflets to assess whether the PVC accessories are compatible and suitable for the intended use.
- The product can be used by multiple people, however it is recommended that every user should use their own set of accessories. In addition please make sure to disinfect the whole set in between use, otherwise it may cause bacterial infection.
- This product should be purchased and used under the guidance of a doctor.
- Do not use highly concentrated medicinal solutions. Please consult your doctor for specific types of nebulization medication, and follow doctor's recommendations.

- Charge the device as soon as low power is indicated.
- If the device will not be used for an extended period of time, please make sure to charge it periodically.
- Ensure that a guardian is present when the nebulizer is used by children under the age of 12.
- Do not store or carry the device with medication / liquid in the medication cup.
- Disposal of waste main units and accessories shall follow the local government regulations.
- The use of this product is different from the laryngeal and nasal mucosa humidification equipment.
- This product must not be used in respiratory anesthesia systems and ventilator systems.
- The service life of this product is 10 years (consumables are excluded).
- Accessories supplied are single-use accessories, which have been sterilized by ethylene oxide. Remember to check the package before use, do not use any accessory supplied in a damaged package but contact the supplier.
- Circuit diagram, list of components as well as more detailed information needed for maintenance is available on request

2. General

2.1 Function and application

Intended use:

The product is intended to aerosolize liquid medications for inhalation by the patient. The device may be used by children and adult patients at home, at the hospital and at

6 sub-acute care settings.

Indication:	Atomization treatment.
Patient population:	The nebulizer can be used under supervision by people over 2 years of age and for self-treatment by people over 12 years of age.
Intended users:	Professional medical staff and patients under their guidance.
Contraindications:	The product must not be used with Pentamidine.

2.2 Specifications

Power supply:	DC 5V or 3.7V rechargeable lithium battery
Input power:	<15 VA
Nebulization rate in level I:	≥ 0.25mL/min
Nebulization rate in level II:	≥ 0.15mL/min
Noise:	≤50 dB
Weight:	0.1kg
Equivalent volume particle diameter distribution:	the volume distribution ratio of small atomized particles (diameter <5 μm) is no less than 60 %
Type of protection against electric shock:	Class II equipment, internally powered equipment
Degree of protection against electric shock:	type BF applied part
Degree of protection against ingress of liquid:	IP22

Note: Please purchase medical power adapter from qualified manufacturer (input: AC 100-240V, 50Hz-60Hz, output: DC 5V, 1.0A).

2.3 Ambient conditions during operation

Temperature:	5 °C - 40 °C
Humidity:	15 % - 90 %
Atmospheric pressure:	700 hPa - 1060 hPa

Attention:

- The device is not suitable for use in a strong electromagnetic interference environment, such as medium and high frequency therapeutic equipment, transformers or large electric cabinet, TV transmission towers, other radio frequency transmission equipment, other electrical appliances that may cause interference, etc.
- When the transport and storage temperature exceed the operating temperature range, the device should be placed in a normal temperature environment for more than 4 hours before use.

2.4 Principles

The Mesh Nebulizer uses a piezoelectric ceramic actuator (PZT) to make a micro-perforated mesh (membrane) vibrate. Whenever a liquid is in contact with the vibrating membrane, it is "extruded" into a cloud of micro-droplets resulting in nebulization and generation of a mist or spray. When the patient uses the mask to cover his/her nose and mouth, the atomized medication is breathed in through his/her inhalation.

Principle of treatment:

The respiratory system is an open system. The atomized medication, after inhalation, can be directly adsorbed on the patient's oral cavity, throat, trachea, bronchus and pulmonary alveoli, etc., through mucous membrane absorption thus achieving the


8 purpose of treatment.

3. Product composition

EN

Structure: The nebulizer set consists of a body, a medication cup, two masks (adult & child), a mouth piece, a power cord, a bag and an adapter (optional).

Main unit:

Accessories:

Adult mask

Child mask

Mouthpiece

Adapter (optional)

Power cord

4. How to use

4.1 Assembly

1. Remove all packages.

Attention:


- Please clean and disinfect the device before using it for the first time.

2. Assembly of the main unit.

1. Fit the medication cup to the main unit by pushing it towards the main unit.

Attention:


- When installing the medication cup to the device body, please make sure to install it to its correct position, indicated by a clicking sound. Incorrect assembly may cause electrode conduction failure thus leading to malfunctioning of the device.

2. In order to remove the medication cup from the main unit, please press and hold the "PUSH" button on the main unit while at the same time pushing the medication cup away from the main unit.

Attention:

- In order to avoid device damage, please press the "PUSH" button first when removing medication cup.

4.2 Instructions for use during treatment

Preparations before use:

1. Remove the medication cup, clean and disinfect it before use.
2. Infusion of medication: Release the clip, lift the cup lid and fill in the medication, as shown below:

Attention:



- Before using any pharmaceutical products or medicines, please consult your doctor to ensure that you are using the product correctly.
- Do not use any highly concentrated, highly viscous, oily or volatile medication. Otherwise this could lead to abnormal atomizing.
- It is recommended not to exceed the maximum capacity of the medication cup. If there is medication in the medication cup, be sure to close the cup lid to prevent leakage.

3. Close the cup lid and fasten the clip.

4. Install the medication cup to the device body.

5. Assemble the mask, as shown below:

Attention:

- If it is necessary to use a new mask for another cycle of inhalation, please repeat the above operation methods to continue using.

4.3 Indicator description

The green indicator (fast breathing) flickers.	Operating in level I
The green indicator (slow breathing) flickers.	Operating in level II
The green indicator glows continuously.	Charge complete
The blue indicator glows continuously.	Charging
The orange indicator flickers.	Low battery
The green indicator flickers at regular intervals.	Cleaning mode

Operation method:

1. **Startup:** short press ON/OFF button to turn on the device, the green light starts flickering and the device starts atomizing.
2. **Mode switch:** the device has two working modes: level I and level II. Level I is the high nebulization rate mode, level II is the low nebulization rate mode. To switch from level I to level II please press the ON/OFF button.
3. **Inhalation:** It is important to use the correct breathing technique to ensure that the particles are distributed as widely as possible across the airways. To ensure that the particles reach your airways and lungs, you must breathe in slowly and deeply, then hold your breath briefly (5 to 10 seconds) and then breathe out quickly.
4. **Shutdown:** Inhaling at level I, please press the ON/OFF button twice to turn off the nebulizer (nebulization stops, indicator goes out). Inhaling at level II you only need to press the ON/OFF button once to turn off the device.

Attention:

1. There is an air hole in the lid of the medication cup, do not cover it to ensure normal atomization.
2. During use, please hold the device steady and do not shake it strongly.
3. After pressing the ON/OFF button, there is a short startup period (within 1 s) before the device actually starts atomizing.
4. When the liquid in the medication cup is nearly used up, it is recommended to slightly tilt the device towards the user (the front side facing the user), so that the remaining solution contacts the nebulizing sheet for ongoing atomization.
5. Due to the different characteristics of medications, the device may not shutdown automatically when some medication is used up. In this case the device needs to be turned off manually to protect the nebulizing mesh from damage.
6. If there is not any medication or liquid left in the medication cup, the device will automatically shut down. Automatic shutdown is indicated by the orange light flickering several times. It is absolutely normal that there is some residual of medication in the medication cup after automatic shutdown.
7. No cycle of inhalation should last any longer than 20 minutes. If you experience any discomfort during use, please turn off the device immediately and consult a doctor.

4.4 Use and Charging of lithium battery**Battery use period:**

1. The indicator lights up in orange as soon as the battery is low, prompting the user to charge the device.
2. Under proper conditions of use, the battery can continuously work for 1 hour after full charge.

Charging:

1. Insert the DC-end of the adapter into the power interface of the device.
2. Plug the adapter into a power socket or use a power bank for charging.

Attention:

- Specification of power adapter: input: AC100-240V 50-60Hz, output: DC5V, 1.0A.
- Please unplug the power adapter after use, avoid staying connected over an extended period of time.
- When the device is left unused for a long time, it should be charged every six months, which will greatly extend the battery life. Batteries cannot be replaced by the user itself. If you need to replace the battery please contact your local service center or our company.
- If you choose to use a different adapter, it should meet the requirements in IEC60601-1, and its input shall be AC100-240V 50-60Hz, output shall be DC5V, 1.0A.

Life of lithium battery:

Under normal conditions, 60 % of the battery capacity should remain after repeated use of about 200 times.

5. Maintenance, Transport and Storage

5.1 Cleaning and disinfection

Clean and disinfect the device after each use. If the device is not cleaned, drying and solidification of the medication will cause aging of both the electrode and the nebulizing mesh, which will affect proper atomizing.

1. Remove the medication cup and any accessory from the device body.
2. Open the cup lid and discard residual medication.
3. Pour 75% ethanol solution into the medication cup, close the cup lid, then leave it for at least 10mins; it is recommended to gently shake it for better disinfection.
4. Immerse the accessories to be disinfected in 75% ethanol solution in a container with lid, close the lid. Let it soak for 10mins or longer.
5. Empty the medication cup and remove the accessories from the disinfectant; clean the medication cup and accessories with clear-water repeatedly.
6. Fill the medication cup with clear-water, fit it to the device body, long press the ON/OFF button (for at least 3 seconds) to enter the cleaning mode and let the device work for about 10mins in order to clean the nebulizing mesh.
7. After cleaning, use medical gauze to absorb any remaining water.
8. Use 75% medicinal alcohol to wipe the surface of device body clean, then air-dry or wipe-dry with a clean, soft cloth.
9. After completion of all cleaning steps, please store the device body, the medication cup and all accessories in a dry, clean place.

Attention:

- Please turn off the device during cleaning and disinfection, do not connect it to a power supply.
- Do not use purified water or distilled water for atomization.
- Please do not maintain the device during normal operation.
- Do not use boiling water for disinfection of neither medication cup nor accessories, otherwise it could lead to deformation of these items.
- All parts disinfected with disinfectant must be meticulously cleaned, or the residual disinfectant may cause symptomatic deterioration.

5.2 Medication cup replacement

The nebulizing mesh needs to be replaced at regular intervals. Generally, the service life of a nebulizing mesh is half a year (if it is used 3 times a day, 20 mins each time.). Its service life depends on the use method, the type of medication, and the extent of cleaning. If the nebulizer fails to atomize or shows only reduced atomizing, please replace the medication cup in time. (If you need to purchase a medication cup, please contact your local retailer.)

5.3 Transport and storage

Ambient conditions during transport and storage:

Temperature:	-40 °C - +55 °C
Relative humidity:	5 % - 96 %
Atmospheric pressure:	500 hPa - 1060 hPa

Transport and storage requirements:

- Do not expose the nebulizer to corrosive gas and keep it in well-ventilated rooms.
- When storing the device, keep it away from children, pets and insects to avoid affecting its performance.
- Do not expose the device to direct sunlight, high temperatures, above average humidity, dust, cotton wool or splash water, etc.
- Avoid violent shock, vibration, rain and snow splash during transportation or storage.

5.4 Pollution-free disposal and recycle

The service life of this product is 10 years.

1. Devices no longer to be used can be returned to the manufacturer or distributor for proper recycling.
2. Used parts can be returned to the manufacturer or distributor for disposal, or they may be disposed of by the user in accordance with relevant laws and regulations.

6. Troubleshooting

Problems	Reason analysis	Solutions
The device does not start.	Low battery.	Please charge the device.
No atomizing or little atomizing during operation.	Medication cup is not well installed.	Check the installation of the medication cup and reinstall it.
	No medication in medication cup.	Pour medication into the medication cup, remember to not exceed its maximum capacity.
	Improper medication.	Consult a doctor to verify that the medication is suitable for the device.
	The nebulizing sheet is clogged.	Clean the medication cup.
There is water around the nozzle of the nebulizer.	Due to temperature differences medication mist might condense when getting in contact with the nozzle.	Remove medication cup from the main body and pour out the water.

After startup, power indicator lights about 1s, then immediately goes out.	Medication cup is not well installed.	Check the installation of the medication cup and reinstall it.
	No medication in medication cup.	Pour medication into the medication cup, remember to not exceed its maximum capacity.
After startup, power indicator lights up, but immediately goes out, or the device does not work properly at all.	Battery is dead.	Please charge the device.
Nebulizer doesn't automatically shut down when medication is used up.	There might be air bubbles in the medication container, which are preventing continuous contact with the mesh.	Press ON/OFF button to turn off the device, and remove the bubbles.
	There might still be some residual medication on the nebulizing mesh.	Press ON/OFF button to turn off the device, and clean the medication cup.
	The electrodes of the medication cup or on the main body may be dirty.	Press ON/OFF button to turn off the device, and clean the electrodes.
The battery use period is too short after the device has been charged.	Battery has not been fully charged.	Please charge the device.
	Battery is damaged.	Please contact your local customer service.
If the device still does not work properly after following this troubleshooting guide, please contact our after-sales service.		

7. Symbols

Symbol	Meaning
IP22	Protected against foreign objects ≥ 12.5 mm and against water dripping at an angle $\leq 15^\circ$.
	Class II equipment
	Type BF applied part
	Keep dry
	Fragile, handle with care
	This way up
	Humidity limitation: 5 % - 96 %
	Temperature limitation: -40°C - $+55^\circ\text{C}$
	Atmospheric limitation: 500 hPa - 1060 hPa
	Refer to instruction manual/booklet Note: On ME equipment "follow instructions for use"
	Serial number
	This item is compliant with Medical Device Directive 93/42/EEC of June 14, 1993, a directive of the European Economic Community.
	Manufacturer

	Do not re-use.
	Consult instructions for use.
	Sterilized using ethylene oxide
	Sterile
	Catalogue number
	Do not use if package is damaged.
	Date of manufacture
	Use-by date
	Batch code
	Authorized representative in the European Community
	Indoor use
	WEEE (2012/19/EU)
	Stacking layer limit is N N subject to actual conditions
	Medical Device

8. Packing List

1x Device body

1x Medication cup

1x Accessories set (adult mask, child mask, mouthpiece)


1x Bag

1x USB cable

1x User manual

Appendix I

Curve chart of equivalent volume particle diameter distribution:

22 The median particle diameter ($D_{0.50}$) is 1-4 µm. Error shall be within $\pm 25\%$.

Appendix II EMC Guidance and Manufacturer Declaration

EN

Warnings:

- Don't near active HF SURGICAL EQUIPMENT and the RF shielded room of an ME SYSTEM for magnetic resonance imaging, where the intensity of EM DISTURBANCES is high.
- Use of the NE-M03C adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, NE-M03C and the other equipment should be observed to verify that they are operating normal.
- Use of accessories, transducers and cables other than those specified or provided by the manufacturer of NE-M03C could result in increased electromagnetic emissions or decreased electromagnetic immunity of NE-M03C and result in improper operation.
- Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the NE-M03C, including cables specified by the manufacturer. Otherwise, degradation of the performance of NE-M03C could result.

Note:

- NE-M03C needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided below.
- The basic performance: Nebulization rate: ≥ 0.25 mL/min.
- When the device is disturbed, the data measured may fluctuate, please measure repeatedly or in another environment to ensure its accuracy.
- Other devices may affect this device even though they meet the requirements of CISPR.

- The following cable types must be used to ensure that they comply with interference radiation and immunity standards:

Name	Cable length (m)
Power cord	1.0

Table 1:

Guidance and Declaration - Electromagnetic Emissions	
Emissions test	Compliance
Radiated RF EMISSIONS CISPR 11	Group 1
Radiated RF EMISSIONS CISPR 11	Class A
Harmonic distortion IEC 61000-3-2	Class A
Voltage fluctuations and flicker IEC 61000-3-3	Complies

Table 2

Guidance and Declaration - Electromagnetic Immunity		
Immunity Test	IEC 60601 Test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	±8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines	±2 kV for power supply lines

Surge IEC 61000-4-5	$\pm 0,5 \text{ kV}, \pm 1 \text{ kV}$, line(s) to line(s)	$\pm 0,5 \text{ kV}, \pm 1 \text{ kV}$, line(s) to line(s)
Voltage dips and voltage interruptions IEC 61000-4-11	0 % UT; 0,5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°. 0 % UT; 1 cycle and 70 % UT; 25/30 cycles; Single phase: at 0°. 0 % UT; 250/300 cycle	0 % UT; 0,5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°. 0 % UT; 1 cycle and 70 % UT; 25/30 cycles; Single phase: at 0°. 0 % UT; 250/300 cycle
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	30 A/m 50Hz/60Hz	30 A/m 50Hz/60Hz
Conducted RF IEC61000-4-6	3 V 0,15MHz - 80 MHz 6 V in ISM bands between 0,15MHz to 80 MHz 80%AM at 1kHz	3 V 0,15MHz - 80 MHz 6 V in ISM bands between 0,15MHz to 80 MHz 80%AM at 1kHz
Radiated RF IEC61000-4-3	10V/m 80 MHz - 2,7GHz 80%AM at 1kHz	10V/m 80 MHz - 2,7GHz 80%AM at 1kHz
NOTE: UT is the a.c. mains voltage prior to application of the test level		

Table 3

Guidance and manufacturer's declaration - electromagnetic Immunity						
Radiated RF IEC61000-4-3 (Test specifications for ENCLOSURE PORT IMMUNITY to RF wireless communications equipment)	Test Fre- quency (MHz)	Band (MHz)	Service	Modula- tion	IEC60601-1-2 Test level (V/m)	Compliance level (V/m)
	385	380 –390	TETRA 400	Pulse modula- tion b) 18 Hz	27	27
	450	430- 470	GMRS 460, FRS 460	FM c) ± 5 kHz deviation 1 kHz sine	28	28
	710	704 – 787	LTE Band 13, 17	Pulse modula- tion b) 217 Hz	9	9
	745					
	780					
	810	800 – 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modula- tion b) 18 Hz	28	28
	870					
	930					

	1720	1700 – 1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation b) 217 Hz	28	28
	1845					
	1970					
	2450	2400 – 2570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation b) 217 Hz	28	28
	5240	5100 – 5800	WLAN 802.11 a/n	Pulse modulation b) 217 Hz	9	9
	5500					
	5785					

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Prolinx GmbH
Contec Medical Systems Co., Ltd.